


CREPE PAPER



Directions For Use

The materials which is to be sterilized are placed into the wire basket.

Wire basket is wrapped by crepe paper.

After sterilization process if any tear or damage on the product please repeat process once more.

Crepe paper is single use only.

It can be used for steam and ethylene oxide, formaldehyde sterilization process.

Interpretation

If necessary, the indicator is placed into the wire basket.

After wrapping process autoclave tape can be used to prevent crepe paper opening.

There are two kind of color. These colors are blue and green.

Storage

Keep products at room temperature 25°C

Expiry Date

The expiry date is printed on the packaging label.

Do not use after expiration date listed on the packs.


Lot Number

Unique identification code, LOT is printed on each packaging labels.

Interfering Substances or Conditions

There are no interfering substances or conditions that could affect the intended use of the crepe paper or adversely affect the crepe paper performance.

CREPE PAPER



Release of Toxic Substances

The crepe paper releases no toxic substances in sufficient quantities to cause either a health hazard or a hazard to the intended properties of the product being sterilized during or after the sterilization process.

Additional Information

Crepe paper is used for wrapping materials.

Steam sterilization process is made at 134°C for 3,5min and 121°C for 15min

Ethylene Oxide sterilization process is made at 40-500 C.

Formaldehyde Sterilization process is made at 50-600 C.

Shelf life of crepe paper is 5 years. The expiry date should be checked before use.

Do not use damaged products.

Crepe Paper Packaging

Product	Weight of Box (kg)	Pcs in a Box	Box Dimension(cm)		
Crepe Paper 50cm x 50cm	8	500	60	55	10
Crepe Paper 60cm x 60cm	12.5	500	65	65	10
Crepe Paper 75cm x 75cm	9.5	250	80	50	11
Crepe Paper 90cm x 90cm	13.5	250	100	53	10
Crepe Paper 100cm x 100cm	17	250	107	60	13
Crepe Paper 120cm x 120cm	12	125	128	52	8