

TYVEK REEL



Indication for Use

Tyvek reels should be cut according to the size of material/device which will be exposed to the sterilization cycle.

The material/device is placed into Tyvek Reel. The open end is sealed between 130-140 degrees with suitable heat sealer.

If necessary, chemical indicator strip can be put into the package.

Color change

After sterilization color change as;

Red to Yellow

Interpretation

Tyvek reels are 5, 7.5, 10, 15, 20, 25, 30, 35 and 40 cm (+/-1 cm) width and 70m. length (+/-1 m)

Storage

Keep products at room temperature 25°C

Expiry Date

Do not use after expiration date listed on the packs.

Lot Number

Unique identification code, LOT is printed on each reel.

Interfering Substances or Conditions

There are no interfering substances or conditions that could affect the intended use of the indicator or adversely affect the reel performance.

Release of Toxic Substances

The indicator releases no toxic substances in sufficient quantities to cause either a health hazard or a hazard to the intended properties of the product being sterilized during or after the sterilization process.

TYVEK REEL



Additional Information

Tyvek Reel will be ready to sterilization cycle after checking the sealing quality.

Please pay attention to the deformation and tear on the packs.

Do not overload the packs.

Please pay attention to the seal strength of both ends.

Do not use damaged products

It is suitable for hydrogen peroxide sterilization.

It has one indicator on it. The indicator is for plasma process.

It is suitable for ethylene and steam sterilization process as well.

But there is no EO and FO indicators on it.

Tyvek Reel Packaging

Dimension	Weight of Box (kg)	Pcs in a Box	Box Dimension(cm)		
5 cm x 70m	5.5	10	42	23	27
7,5 cm x 70m	6.5	8	42	23	32
10 cm x 70m	8.2	8	42	23	42
15 cm x 70m	6.5	4	42	23	32
20 cm x 70m	8	4	42	23	42
25 cm x 70m	5.2	2	42	23	27
30 cm x 70m	6	2	42	23	32
35 cm x 70m	7.5	2	42	23	37
40 cm x 70m	8.5	2	42	23	42