

SECURITY SEAL



Directions For Use

4A Security Seals are used to secure the container for opening by un-authorized staff. Container can be opened only by breaking the seal. Security seals are guaranteeing the sterility of materials from the moment of sterilization until the actual use. 4A is offering Security Seals with and without indicators.

Color change

If colors of all lines change from blue to green the operation process is considered successfully.

Safety Precautions

Do not use damaged products.

Storage

Keep products at room temperature 25°C

Expiry Date

The expiry date is printed on the label. Do not use after expiration date listed on the packs.

Lot Number

Unique identification code, LOT is printed on each packaging labels.

Interfering Substances or Conditions

There are no interfering substances or conditions that could affect the intended use of the indicator or adversely affect the indicator performance.

Release of Toxic Substances

The indicator releases no toxic substances in sufficient quantities to cause either a health hazard or a hazard to the intended properties of the product being sterilized during or after the sterilization process.